


YOU CAN DO IT WITHOUT FINGERSTICKS*

The FreeStyle Libre 14 day system is an accurate,¹ convenient, and easy² way to get glucose readings and trends.



Two components:

innovative sensor technology




- Start taking glucose readings just one hour after sensor application*
- Small, discreet,² and water resistant†
- Can be worn for up to 14 days
- No fingerstick calibration required
- Sensor automatically measures glucose readings day and night‡

compact, lightweight reader



- Painless,² one-second scan—even through clothing§
- Stores up to 90 days of glucose data
- Backlit touchscreen for visibility



How to get started with the FreeStyle Libre 14 day system

1 get a prescription
Request a written prescription. Your healthcare professional can use our example as a guide.

2 go to pharmacy
Fill your prescription at a participating pharmacy.[¶] Your pharmacy may have to order the product for you.

FreeStyle Libre 14 day system is available with a prescription at participating pharmacies.^{¶¶}

FREESTYLE LIBRE INDICATIONS AND IMPORTANT SAFETY INFORMATION

The FreeStyle Libre 14 day Flash Glucose Monitoring System is a continuous glucose monitoring (CGM) device indicated for the management of diabetes in persons age 18 and older. It is designed to replace blood glucose testing for diabetes treatment decisions. The System detects trends and tracks patterns aiding in the detection of episodes of hyperglycemia and hypoglycemia, facilitating both acute and long-term therapy adjustments. Interpretation of the System readings should be based on the glucose trends and several sequential readings over time. The System is intended for single patient use and requires a prescription.

CONTRAINDICATIONS: Remove the sensor before MRI, CT scan, X-ray, or diathermy treatment.

WARNINGS/LIMITATIONS: Do not ignore symptoms that may be due to low or high blood glucose, hypoglycemic unawareness, or dehydration. Check sensor glucose readings with a blood glucose meter when Check Blood Glucose symbol appears, when symptoms do not match system readings, or when readings are suspected to be inaccurate. The FreeStyle Libre 14 day system does not have alarms unless the sensor is scanned, and the system contains small parts that may be dangerous if swallowed. The FreeStyle Libre 14 day system is not approved for pregnant women, persons on dialysis, or critically-ill population. Sensor placement is not approved for sites other than the back of the arm and standard precautions for transmission of blood borne pathogens should be taken. The built-in blood glucose meter is not for use on dehydrated, hypotensive, in shock, hyperglycemic-hyperosmolar state, with or without ketosis, neonates, critically-ill patients, or for diagnosis or screening of diabetes. Review all product information before use or contact Abbott toll-free 855-632-8658 or visit www.FreeStyleLibre.us for detailed indications for use and safety information. For full indications for use and safety information, visit <https://www.FreeStyleLibre.us/safety-information.html>.

The product images are for illustrative purposes only.

* Fingersticks are required for treatment decisions when you see Check Blood Glucose symbol, when symptoms do not match system readings, when you suspect readings may be inaccurate, or when you experience symptoms that may be due to high or low blood glucose.

† Sensor is water resistant in up to 1 meter (3 feet) of water. Do not immerse longer than 30 minutes.

‡ Based on the sensor being replaced once every 14 days, and scanned at least once every 8 hours.

§ The reader can capture data from the sensor when it is within 1 cm to 4 cm of the sensor.

¶ For Medicare patients, please visit FreeStyleLibre.us/buying-guide for more information.

¶ Participating pharmacies are subject to change without notice.

References: 1. FreeStyle Libre User Manual. 2. Data on File: Abbott Diabetes Care.

FreeStyle, Libre, and related brand marks are trademarks of Abbott Diabetes Care Inc. in various jurisdictions. Other trademarks are the property of their respective owners.

©2018 Abbott. ADC-05571 v7.0 10/2018